

DEUTERATED NMR SOLVENTS-HANDY REFERENCE DATA

Compound Mol. Wt.	d_4^{20}	m.p.*	b.p.*	$d_w(\text{mult})^*$	J_{HD}	$d_c(\text{mult})^*$	$J_{CD}(J_{CF})$
Acetic Acid-d ₄ 64.078	1.12	17	118	11.53 (1) 2.03 (5)	2	178.4 (br) 20.0 (7)	20
Acetone-d ₆ 64.117	0.87	-94	57	2.04 (5)	2.2	206.0 (13) 29.8 (7)	0.9 20
Acetonitrile-d ₃ 44.071	0.84	-45	82	1.93 (5)	2.5	118.2 (br) 1.3 (7)	21
Benzene-d ₆ 84.152	0.95	5	80	7.15 (br)		128.0 (3)	24
Chloroform-d 120.384	1.50	-64	62	7.24 (1)		77.0 (3)	32
Cyclohexane-d ₁₂ 96.236	0.89	6	81	1.38 (br)		26.4 (5)	19
Deuterium Oxide 20.028	1.11	3.8	101.4	4.63(DSS) 4.67(TSP)			
1,2-Dichloroethane-d ₂ 102.985	1.25	-40	84	3.72 (br)		43.6 (5)	23.5
Diethyl-d ₁₀ Ether 84.185	0.82	-116	35	3.34 (m) 1.07 (m)		65.3 (5) 14.5 (7)	21 19
Diglyme-d ₁₂ 148.263	0.95	-68	162	3.49 (br) 3.40 (br) 3.22 (5)	1.5	70.7 (5) 70.0 (5) 57.7 (7)	21 21 21
Dimethylformamide-d ₂ 80.138	1.04	-61	153	8.01 (br) 2.91 (5) 2.74 (5)	2 2	162.7 (3) 35.2 (7) 30.1 (7)	30 21 21
Dimethyl-d ₄ Sulphoxide 84.170	1.18	18	189	2.49 (5)	1.7	39.5 (7)	21
p-Dioxane-d ₈ 96.156	1.13	12	101	3.53 (m)		66.5 (5)	22
Ethyl Alcohol-d ₆ (anh) 52.106	0.91	<-130	79	5.19 (1) 3.55 (br) 1.11 (m)		56.8 (5) 17.2 (7)	22 19
Glyme-d ₁₀ 100.184	0.86	-58	83	3.40 (m) 3.22 (5)	1.6	71.7 (5) 57.8 (7)	21 21
Hexafluoroacetone Deuterate* 198.067	1.71	21		5.26 (1)		122.5 (4) 92.9 (7)	(287) (34.5)
HMPT-d ₁₂ 197.314	1.14	7	106 (11)	2.53(2 x 5)	2(9.5)	35.8 (7)	21
Methyl Alcohol-d ₄ 36.067	0.89	-98	65	4.78 (1) 3.30 (5)	1.7	49.0 (7)	21.5
Methylene Chloride-d ₂ 86.945	1.35	-95	40	5.32 (3)	1	53.8 (5)	27
Nitrobenzene-d ₅ 128.143	1.25	6	211	8.11 (br) 7.67 (br) 7.50 (br)		148.6 (1) 134.8 (3) 129.5 (3) 123.5 (3)	24.5 (p) 25 26
Nitromethane-d ₃ 64.059	1.20	-29	101	4.33 (5)	2	62.8 (7)	22
isoPropyl Alcohol-d ₈ 68.146	0.90	-86	83	5.12 (1) 3.89 (br) 1.10 (br)		62.9 (3) 24.2 (7)	21.5 19
Pyridine-d ₅ 84.133	1.05	-42	116	8.71 (br) 7.55 (br) 7.19 (br)		149.9 (3) 135.5 (3) 123.5 (3)	27.5 24.5 (y) 25
Tetrahydrofuran-d ₂ 80.157	0.99	-109	66	3.58 (br) 1.73 (br)		67.4 (5) 25.3 (br)	22 20.5
Toluene-d ₈ 100.191	0.94	-95	111	7.09 (m) 7.00 (br) 6.98 (m) 2.09 (5)	2.3	137.5 (1) 128.9 (3) 128.0 (3) 125.2 (3) 20.4 (7)	23 24 24 (p) 19
Trifluoroacetic Acid-d ⁴ 115.030	1.50	-15	72	11.50 (1)		164.2 (4) 116.6 (4)	(44) (283)
2,2,2-Trifluoroethyl Alcohol-d ⁴ 103.059	1.45	-44	75	5.02 (1) 3.88(4 x 3)	2 (9)	126.3 (4) 61.5(4 x 5)	(277) 22 (36)

*Melting and boiling points (in °C) are those of the corresponding light compound (except for D₂O) and are intended only to indicate the useful liquid range of the materials.

¹H (of the residual protons) and ¹³C spectra were determined on HA-100 and XL-100-15 spectrometers, respectively, for the same sample of each solvent containing 5% TMS (v/v). The chemical shifts are in ppm relative to TMS; the coupling constants are in Hz. (Since deuterium has a spin of 1, triplets arising from coupling to deuterium have the intensity ratio of 1:1:1, etc.) The multiplicity br indicates a broad peak without resolvable fine structure, while m denotes one with fine structure. It should be noted that the chemical shifts, in particular, can be dependent on solute, concentration and temperature.

P.O. BOX 2951 TERMINAL ANNEX
LOS ANGELES, CA 90051
213 723-9521

800 325-9034

4545 OLEATHA AVE.
ST. LOUIS, MO 63116
314 353-7000

MERCK & CO., Inc.



ISOTOPES

P.O. BOX 899
POINTE CLAIRE-DORVAL
QUEBEC H9R 4P7
514 697-2823

MERCK
SHARP
& DOHME CANADA LIMITED
ISOTOPES



**MERCK
SHARP
& DOHME** CANADA LIMITED
ISOTOPE DIVISION

ISOTOPICS

P.O. Box 899, Pointe Claire-Dorval H9R 4P7, Quebec - Tel. 514-697-2823 - Cable: ISOMERCK/Telex: 05-821-533/TWX: 610-421-3617

SUGGESTIONS FOR USE & HANDLING OF NMR SOLVENTS

In Fourier Transform NMR studies, especially when using microgram quantities, peaks arising from water and hydrocarbons leached from cap liners may present serious problems. These extraneous peaks may be avoided by the use of solvents sealed in glass ampoules.

Our NMR solvents are carefully dried prior to packaging. When opened, most solvents will absorb atmospheric moisture. To minimize the water-absorption problem we recommend the use of our sealed glass ampoules, available in 1 gram or 0.5 ml sizes. They are especially convenient for glove-box use. For a number of solvents there is no extra charge for 10 gram quantities - i.e., the price for 10 x 1 gram is the same as for 1 x 10 grams.

We find that solvents can be rapidly dried with Linde 3A or 4A sieves, 1/16"; bead-shaped sieves do not perform satisfactorily, even when freshly conditioned.

To avoid the introduction of water by desorption from the wall of the NMR tube (especially important when overnight runs are required), the tubes can be preconditioned with D_2O . This does not remove the water, but substitutes D_2O for the signal-producing H_2O which, for example, gives rise to a peak at ca. δ 1.5 in $CDCl_3$. An easy method is to keep one ml or more of D_2O in the NMR tube overnight, decant, wash with a small amount of CD_3OD or $(CD_3)_2CO$ to remove the residual D_2O , dry briefly with a stream of dry nitrogen and store in a desiccator saturated with D_2O . Using this technique, a freshly-dried (Linde 3A sieve, 1/16") $CDCl_3$ sample analysed in a D_2O conditioned tube leaves the water peak sufficiently small for most purposes.

Our laboratories have tested the use of septum vials for preserving "100%" D_2O and have found that with a little additional care, decrease in isotopic enrichment can be avoided even with the repeated use of a syringe to withdraw material. We suggest that first the syringe should be flushed with dry nitrogen; then a quantity of dry nitrogen equivalent to the amount of D_2O required should be injected into the septum vial before the material is withdrawn. In our stability tests, we have observed no decrease in isotopic enrichment.

If the cap is kept properly tightened, TMS will not readily evaporate at room temperature. If kept refrigerated, ensure that the TMS bottle is at room temperature before it is opened in order to avoid condensation of atmospheric water in the cold bottle. If two layers are present in the TMS bottle, the lower one will be the water.

Proton Chemical Shifts
of Common Impurities

Impurity	Solvent				
	CDCl ₃	(CD ₃) ₂ SO	C ₅ D ₅ N	C ₆ D ₆	D ₂ O
Acetic Acid	2.13	1.95	2.13	1.63	2.16
Acetone	2.17	2.12	2.00	1.62	2.22
Acetonitrile	1.98	2.09	1.85	0.67	2.05
Amisyl Alcohol	4.61 3.79	4.44 3.74	4.92 3.68	4.36 3.32	4.57 3.81
Benzene	7.37	7.40	7.33	7.30	7.44
Bromoform	6.85	7.75	7.86	5.89	*
n-Butanol	3.67(t) 0.94(t)	3.41(t) 0.89(t)	3.80(t) 0.89(t)	3.38(t) 0.83(t)	3.60 0.89(t)
t-Butanol	1.28	1.14	1.37	1.06	1.23
Chloroacetic Acid	4.14	4.28	4.44	3.31	4.25
Chloroform	7.27	8.35	8.41	6.41	*
Cyclohexane	1.43	1.42	1.38	1.40	*
1,2 Dibromoethane	3.63	3.84	3.70	2.88	3.79
Diochloroacetic Acid	5.98	6.68	6.72	5.43	6.21
1,2 Dichloroethane	3.73	3.93	3.78	2.99	3.92
Dichloromethane	5.30	5.79	5.62	4.46	*
Diethyleneglycol dimethyl ether	3.60(m) 3.38	3.49(m) 3.28	3.57(m) 3.27	3.40(m) 3.13	3.63(m) 3.37

Shifts are PPM (δ) with respect to TMS.

* Indicates that the impurity is insoluble in the solvent.

+ Indicates that the impurity and solvent react.

Impurity	CDCl ₃	(CD ₃) ₂ SO	C ₅ D ₅ N	C ₆ D ₆	D ₂ O
Diethyl ether	3.48(q) 1.2(t)	3.42(q) 1.13(t)	3.38(q) 1.12(t)	3.27(q) 1.10(t)	3.56(q) 1.17(t)
Di-isopropyl ether	1.12(d)	1.04(d)	1.09	0.50(d)	1.12(d)
Dimethylacetamide	2.08 3.01 2.94	2.99 2.82 1.99	2.82 2.70 1.96	2.59 2.11 1.64	3.05 2.89 2.08
Dimethylformamide	8.01 2.95 2.88	7.98 2.92 2.76	2.72 2.66	2.40 1.98	7.91 3.00 2.86
Dimethyl Sulfoxide	2.62	2.52	2.49	1.91	2.70
Dioxan	3.70	3.61	3.61	3.38	3.75
Ethenediol	3.76	3.42	4.01	4.21	3.66
Ethanol	3.72(q) 1.24(t)	3.49(q) 1.09(t)	3.86(q) 1.29(t)	3.39(q) 0.97(t)	3.64(q) 1.16(t)
Ethyl Acetate	4.12(q) 2.04 1.25(t)	4.08(q) 2.02 1.21(t)	4.06(q) 1.94 1.10(t)	3.91(q) 1.68 0.94(t)	4.14(q) 2.08 1.23(t)
Ethyl Formate	8.04 4.22(q) 1.29(t)	8.23 4.17(q) 1.24(t)	8.22 4.14(q) 1.10(t)	7.60 3.83(q) 0.85(t)	8.16 4.28(q) 1.29(t)
Formic Acid	8.02	8.18	8.54	7.24	8.22
Furfuraldehyde	9.68 7.71(m) 7.25(m) 6.60(m)	9.69 8.14(m) 7.59(m) 6.82(m)	9.78 7.87(m) 7.31(m) 6.56(m)	9.36 6.89(m) 6.57(m) 5.82(m)	9.52 7.94(m) 7.59(m) 6.78(m)
Isobutyl methyl ketone	2.12 0.92(d)	2.08 0.88(d)	2.00 0.83(d)	1.61 0.74(d)	2.19 0.89(d)

Shifts are PPM (δ) with respect to TMS.

* Indicates that the impurity is insoluble in the solvent.

+ Indicates that the impurity and solvent react.

Impurity	CDCl ₃	(CD ₃) ₂ SO	C ₅ D ₅ N	C ₆ D ₆	D ₂ O
Isopropyl Acetate	2.02 1.22(d)	2.00 1.21(d)	1.95 1.13(d)	1.69 1.02(d)	*
Isopropyl Alcohol	4.03(m) 1.20(d)	1.06(d)	4.16(m) 1.29(d)	3.76(m) 1.01(d)	1.18(d)
Methanol	3.48	3.20	3.57	3.09	3.35
Methyl Acetate	3.67 2.05	3.61 2.02	3.55 1.92	3.28 1.63	3.68 2.09
Methyl Iodide	2.16	2.21	+	1.47	*
Morphine	3.69(d) 2.85(d)	3.52(d) 2.68(d)	2.67(d) 2.81(d)	3.50(d) 2.49(d)	3.70(d) 2.79(d)
Nitromethane	4.32	4.44	4.39	3.09	4.41
Petroleum spirit	1.28 0.90(t)	1.28 0.89(t)	1.20 0.86(t)	1.22 0.89(t)	*
Potassium Acetate	*	1.60	2.11	*	1.91
Propane-1,2 diol	1.17(d)	1.04(d)	1.36(d)	1.10(d)	1.13(d)
Propanol	3.60(t) 1.60(m) 0.93(t)	1.45(m) 0.87(t)	3.75(t) 1.70(m) 0.97(t)	3.76(t) 1.40(m) 0.80(t)	3.61(t) 1.57(m) 0.89(t)
Propionic Acid	2.42(q) 1.18(t)	2.26(q) 1.03(t)	2.46(q) 1.20(t)	2.02(q) 0.89(t)	2.47(q) 1.10(t)
Pyridine	8.60(m) 7.69(m) 7.28(m)	8.61(m) 7.83(m) 7.40(m)	8.71(m) 7.58(m) 7.21(m)	8.50(m) 7.05(m) 6.70(m)	8.50(m) 7.90(m) 7.47(m)
Succinimide	2.75	2.63	2.64	1.50	2.78
1,1,2,2 Tetrachloroethane	5.96	6.92	6.90	4.96	*

Shifts are PPM (δ) with respect to TMS.

* Indicates that the impurity is insoluble in the solvent.

+ Indicates that the impurity and solvent react.

<u>Impurity</u>	<u>CDCl₃</u>	<u>(CD₃)₂SO</u>	<u>C₅D₅N</u>	<u>C₆D₆</u>	<u>D₂O</u>
THF	3.74(m) 1.85(m)	3.63(m) 1.78(m)	3.67(m) 1.64(m)	3.01(m) 0.87(m)	3.75(m) 1.88(m)
Toluene	7.19 2.34	7.22 2.32	7.22 2.22	7.10 2.13	*
1,1,1 Trichloroethane	1.72	2.80	2.61	1.58	*
Triethylamine	2.56(q) 1.03(t)	2.47(q) 0.99(t)	2.43(q) 0.96(t)	2.40(q) 0.95(t)	2.59(q) 1.02(t)
Trimethyl borate	3.48	3.46	3.58	3.13	+
Trimethyl Phosphate	3.80(d)	3.72(d)	3.71(d)	3.37(d)	3.82(d)

Shifts are PPM (δ) with respect to TMS.

* Indicates that the impurity is insoluble in the solvent.

+ Indicates that the impurity and solvent react.

ACETONE CIL LOT # F - 5153A



CAMB.026
DATE 24-2-88


SF 300.135
SFO 300.130
SY 100.0
Q1 6600.000
SI 16384
TD 16384
SW 4000.000
HZ/PT .488

PW 9.0
RD 8.000
AQ 2.048
RG 10
NS 16

DR 12
FW 5000
DP 12H PO

LB 0.0
HZ/CM 150.059
PPM/CM 500
SR 4929.13

Acetone-d₆
99.9 atom % D
DLM-9



ACETONE CIL LOT # F - 49398



CAMB.029
DATE 24-2-88


SF 300.135
SFO 300.130
SY 100.0
Q1 6600.000
SI 16384
TD 16384
SW 4000.000
HZ/PT .488

PW 9.0
RD 8.000
AQ 2.048
RG 20
NS 16

DR 12
FW 5000
DP 12H PO

LB 0.0
HZ/CM 150.059
PPM/CM 500
SR 4929.13

Acetone-d₆ "100%"
99.96 atom % D
DLM-38



Matching the quality of the best competition is CIL's minimum level of acceptable performance!

Acetonitrile-d₃
99.8 atom % D
DLM-21




CAMB.073
DATE 29-2-88

SF 300.135
SF0 300.130
SY 100.0
O1 6800.000
SI 16384
TD 16384
SW 5000.000
HZ/PT .610

PX 9.0
RD 8.000
AQ 1.638
RG 4
NS 16

DR 12
FW 6300
DP 12H P0

LB 0.0
HZ/CM 150.049
PPM/CM .500
SR 4965.44



Benzene-d₆
99.6 atom % D
DLM-1




CAMB.062
DATE 25-2-88

SF 300.133
SF0 300.130
SY 100.0
O1 5000.000
SI 16384
TD 16384
SW 5000.000
HZ/PT .610

PX 9.0
RD 8.000
AQ 1.638
RG 4
NS 16

DR 12
FW 6300
DP 12H P0

LB 0.0
HZ/CM 150.049
PPM/CM .500
SR 3385.05





CAMB.005
DATE 24-2-88



SF 300.133
SFO 300.130
SY 100.0
O1 5000.000
SI 16384
TD 16384
SW 5000.000
HZ/PT .610

PW 9.0
RD 8.000
AQ 1.638
R6 10
NS 16

DR 12
FW 6300
DP 12H P0

LB 0.0
HZ/CM 150.073
PPM/CM .500
SR 3364.60

Chloroform-d
99.8 atom % D
DLM-7



CAMB.002
DATE 24-2-88



SF 300.133
SFO 300.130
SY 100.0
O1 5000.000
SI 16384
TD 16384
SW 5000.000
HZ/PT .610

PW 9.0
RD 8.000
AQ 1.638
R6 10
NS 16

DR 12
FW 6300
DP 12H P0

LB 0.0
HZ/CM 150.073
PPM/CM .500
SR 3364.60

Chloroform-d "100%"
99.96 atom % D
DLM-29





CAMB.020
DATE 24-2-88


SF 300.135
SFO 300.130
SY 100.0
Q1 6600.000
SI 16384
TD 16384
SW 4000.000
HZ/PT .488

PW 9.0
RD 8.000
AQ 2.048
RG 10
NS 16

DR 12
FW 5000
OP 12H PO

LB 0.0
HZ/CM 150.039
PPM/CM .500
SR 4791.63

Dimethyl-d₆ sulfoxide
99.9 atom % D
DLM-10



Matching the quality of the best competition is CIL's minimum level of acceptable performance!



CAMB.016
DATE 24-2-88


SF 300.135
SFO 300.130
SY 100.0
Q1 6600.000
SI 16384
TD 16384
SW 4000.000
HZ/PT .488

PW 9.0
RD 8.000
AQ 2.048
RG 10
NS 16

DR 12
FW 5000
DP 12H PO

LB 0.0
HZ/CM 150.039
PPM/CM .500
SR 4791.63

Dimethyl-d₆ sulfoxide "100%"
99.96 atom % D
DLM-34



D20 CIL LOT # F - 5624C



CAMB.057
DATE 25-2-88


SF 300.134
SFO 300.130
SY 100.0
O1 6600.000
SI 16384
TD 16384
SM 10000.000
HZ/PT 1.221

PW 9.0
RD 8.000
AQ .819
RG 4
NS 16

DR 12
FM 12500
DP 20H P0

LB 0.0
HZ/CM 150.049
PPM/CM .500
SR 4132.70

Deuterium oxide
99.9 atom % D
DLM-4



DIOXANE CIL LOT # P - 976



CAMB.093
DATE 7-3-88


SF 300.134
SFO 300.130
SY 100.0
O1 6200.000
SI 16384
TD 16384
SM 4000.000
HZ/PT .488

PW 9.0
RD 8.000
AQ 2.048
RG 10
NS 16

DR 12
FM 5000
DP 12H P0

LB 0.0
HZ/CM 150.059
PPM/CM .500
SR 4486.23

p-Dioxane-d8
99 atom % D
DLM-28



Matching the quality of the best competition is CIL's minimum level of acceptable performance!

METHANOL CIL LOT # F - 5057B




CAMB.037
DATE 24-2-88
SF 300.135
SFO 300.130
SY 100.0
O1 6200.000
S1 16384
T0 16384
SW 4000.000
HZ/PT .488

PW 9.0
RD 8.000
AG 2.048
HG 10
NS 1E

OR 12
FW 5000
DP 12H PC

LB 0.0
HZ, CM 150.059
PPM/CM .500
SR 4561.35

Methanol-d₄
99.8 atom % D
DLM-24



METHANOL CIL LOT # F - 5239C




CAMB.042
DATE 25-2-88
SF 300.135
SFO 300.130
SY 100.0
O1 6200.000
S1 16384
T0 16384
SW 4000.000
HZ/PT .488

PW 9.0
RD 8.000
AG 2.048
HG 10
NS 1E

OR 12
FW 5000
DP 12H PC

LB 0.0
HZ, CM 150.059
PPM/CM .500
SR 4561.35

Methanol-d₄ "100%"
99.96 atom % D
DLM-51



Matching the quality of the best competition is CIL's minimum level of acceptable performance!

METHYLENE CHLORIDE CIL LOT # F - 5669B




CAMB.091
DATE 3-3-88
SF 300.134
SFO 300.130
SY 100.0
O1 5800.000
SI 16384
TD 16384
SN 4000.000
HZ/PT .488

PW 9.0
RD 8.000
AQ 2.048
RG 10
NS 16
DR 12
FW 5000
DP 12H P0

LB 0.0
HZ/CM 150.059
PPM/CM .500
SR 3948.00

Methylene Chloride-d₂
99.9 atom % D
DLM-23



PYRIDINE CIL LOT # F - 4991C




CAMB.080
DATE 29-2-88
SF 300.133
SFO 300.130
SY 100.0
O1 5000.000
SI 16384
TD 16384
SN 5000.000
HZ/PT .610

PW 9.0
RD 8.000
AQ 1.638
RG 4
NS 16
DR 12
FW 6300
DP 12H P0

LB 0.0
HZ/CM 150.049
PPM/CM .500
SR 2917.46

Pyridine-d₅
99.5 atom % D
DLM-13



Matching the quality of the best competition is CIL's minimum level of acceptable performance!

TOLUENE CIL LOT # F - 5451A

~~BRUKER~~

CAMB.063
DATE 29-2-88

SF 300.133
SFO 300.130
SY 100.0
O1 5000.000
SI 16384
TD 16384
SN 5000.000
HZ/PT .610

PW 9.0
RD 8.000
AQ 1.638
RG 4
NS 16

DR 12
FX 6300
DP 12H PO

LB 0.0
HZ/CM 150.049
PPM/CM .500
SR 3400.43

Toluene-d₈
99.6 atom % D
DLM-5

